

# TIARIS™ Pfu DNA Polymerase Kit

## Ordering Info

TBK0085, 12.5 U (2.5U/ μL)(sample)

TBK0086, 125 U (2.5U/ μL)

TBK0087, 500 U (2.5U/ μL)

## Description

TIARIS™ Pfu DNA Polymerase Kit is an ideal kit when high fidelity is an important requirement in your PCR. It includes the proofreading enzyme TIARIS™ Pfu DNA polymerase and the highly purity High-Q™ dNTPs. TIARIS™ Pfu DNA polymerase (90 kDa) is a high fidelity polymerase from *Pyrococcus furiosus* produced in *Escherichia coli*. The enzyme has a strong 5'→ 3' polymerase activity and 3'→5' exonuclease (proofreading) activity which corrects nucleotide incorporation errors, thereby increasing fidelity and accuracy of DNA polymerization.

## Features

- High accuracy and fidelity PCR with an error rate of  $1.3 \times 10^{-6}$  errors/bp per cycle<sup>1</sup>.
- Elongation velocity is 0.2~0.4 kb/min (70~75°C).
- Suitable for blunt cloning purposes, the enzyme produces blunt-end PCR fragments.
- Highly thermostable enzyme, which retains 94-99% of its activity after 1 hour at 95°C.
- PCR amplification > 2 kb require optimization.

## Applications

- Broad range of techniques that require high fidelity PCR products such as: cloning, expression of gene of interest, mutagenesis.
- Generation of PCR fragments for blunt cloning.
- Filling 5' protruding ends.

## Kit Components

Components	TBK0086	TBK0087
TIARIS™ Pfu DNA polymerase (2.5 U/ μL)	50 μL	200 μL
Pfu PCR Buffer (10x)	1.5 mL	1.5 mL
MgCl <sub>2</sub> 25 mM	1.5 mL	1.5 mL
High-Q™ dNTPs 10mM TOTAL	0.2mL	1 mL

**Order Info Kit Components:** TIARIS™ Pfu DNA polymerase (TBZ0212) | Pfu PCR Buffer 10x (TBB0314) | MgCl<sub>2</sub> 25 mM (TBR0215) | High-Q™ dNTPs 10 mM TOTAL (TBR0209).

## Storage

Store at -20°C. Shipped in blue ice.

## Unit Definition

One unit of TIARIS™ Pfu DNA polymerase catalyzes the incorporation of 10 nmol of deoxyribonucleotides into a polynucleotide fraction in 30 minutes at 72°C.

## Quality Control

Functionally tested in a 0.5kb PCR amplification (GC 52%).

## Material required (not supplied)

- PCR Grade Water (TBB0303)
- PCR Tubes
- Primers

<sup>1</sup>Cline J, Braman JC, Hogref HH (1996) *Nucleic Acids Res* 24 (8): 3546–3551.

## PROTOCOL

1. Thawing all components on ice. Vortex them and centrifugate.
2. On ice, prepare a mix of the following components, considering the number of samples plus two extra reactions.

Reaction Components	Final Concentration	Volume	Volume
TIARIS™ Pfu PCR Buffer 10x	1 x	2 µL	5 µL
MgCl <sub>2</sub> 25 mM	2.5 mM	2 µL	5 µL
dNTPs Mix 10 mM TOTAL	0.8 mM	1.6 µL	4 µL
Forward Primer (15 pmol/ µL)	0.2-0.75 µM	0.3-1 µL	1.7-2.5 µL
Reverse Primer (15 pmol/ µL)	0.2-0.75 µM	0.3-1 µL	1.7-2.5 µL
TIARIS™ Pfu DNA polymerase (2.5 U/µL)	0.05 U/ µL	0.4 µL	1 µL
Water, molecular biology grade		up 20 µL*	up 50 µL*
DNA template (add in step 4)		*	*
<b>Final Volume</b>		<b>20 µL</b>	<b>50 µL</b>

\* consider volume of template to be added in step 4.

3. Distribute the mix prepared in each PCR tube or well.
4. Add in each tube the DNA sample (plasmid DNA, cDNA: 10-75 ng; gDNA: 100-500 ng). Mix well.
5. Set up thermocycler with the following suggested parameters,

Process	Cycles	Temperature	Time
Denaturation		94-98 °C*	0:40
Annealing	25 - 30 x	T <sub>m</sub> - 5°C	0:40
Extension		72 °C	2:00 per kb
Final Extension	1 x	72 °C	10:00
Conservation	1 x	4 °C	∞

\* Use 98°C, if GC content is high.

6. Store the PCR samples at -20°C.