Q-PLUS™-Probe Real Time qPCR Master Mix (2x)

Ordering Info

TBK0070, 5 reactions (sample)	TBK0072, 400 reactions	TBK0074, 1000 reactions
TBK0071, 200 reactions	TBK0073, 500 reactions	

Description

Q-PLUS™-Probe Real-Time qPCR Master Mix (2X) is an optimized master mix designed for use in real-time quantitative PCR. The kit includes a qPCR master mix presented as a 2× reaction mixture. This mix incorporates all essential components for real-time PCR, including dNTPs, stabilizers, and enhancers, designed for the efficient amplification and detection of DNA in qPCR based on a wide range of probe-based technologies, including Taqman®, Molecular Beacons® and Scorpion probes®. This kit also includes a separate vial of ROX that can be optionally added to the qPCR reaction Mix. The final concentration of ROX will vary depending on each real-time cycler manufacturer's specification.

Features

- Ready to use.
- Enhanced specificity and sensibility, amplifies low copy number targets with reduced non-specific.
- Compatible with fast and standard PCR program.
- High Efficiency in multiplex reactions.

ROX Reference Dye

The passive reference dye ROX is necessary for certain real-time PCR machines as it compensates for non-PCR-related variations in fluorescence detection. The fluorescence emitted from the ROX dye remains constant throughout the real-time PCR process, providing a stable baseline against which PCR-related fluorescent signals can be normalized. As a result, the ROX dye can compensate any differences in fluorescence detection between wells that may arise from slight variations in reaction volume or differences in well position.

Depending on your equipment, prior to use for the first time, add 18 µL ("High ROX") of the ROX reference dye to the 1 mL Q-PLUS™-Probe Real Time qPCR Master Mix, (2x) and vortex briefly. Once ROX has been added, the Master Mix can be used directly or stored at -20°C for up to 1 year.

Kit Components

Components	TBK0071	TBK0072	TBK0073	TBK0074
Q-PLUS™-Probe Real Time qPCR Master Mix (2x)	2 x 1 mL	4 x 1 mL	5 x 1 mL	10 x 1 mL
PCR Grade Water	2 mL	2 x 2 mL	3 x 2 mL	5 x 2 mL
ROX Reference Dye	1 vial	1 vial	1 vial	2 vial

Order Info Kit Components: ROX Reference Dye (TBR0278) | PCR Grade Water, nuclease free (TBB0303).

Storage

Shipped on blue ice. Upon receipt, kit components should be immediately stored at -20°C. Avoid repeated freezing and thawing. Maintain cold when thawed.

Applications

- Real time quantitative PCR to quantify DNA/ cDNA: gene copy number determination, microbial detection, gene expression, etc.
- Genotyping.
- Multiplex or singleplex PCR

Technical Assistance

Please refer any technical questions to support@tiarisbiosciences.com

PROTOCOL

1. Gently vortex and briefly centrifuge Q-PLUS™-Probe Real Time qPCR Master Mix (2x) after thawing.
2. Place a tube on ice and add the following components for each 20 µL reaction. Prepare sufficient master mix for the number of reactions. Consider one or two extras.

Components	Volume	Final Conc.
Q-PLUS™-Probe Real Time qPCR Mastermix (2x)*	10 µL	1x
Forward primer 15 µM (15 pmol/µL)	<2 µL	50-400 nM
Reverse primer 15 µM (15 pmol/µL)	<2 µL	50-400 nM
Probe (10µM)	0.2-0.8 µL	100-400 nM
Template DNA**(step 5)	<2 µL	**
PCR Grade Water, nuclease-free	up 20 µL	
Final Volume	20 µL	

* ROX as required; see section: "ROX Reference Dye"

**Optimal amounts of template: genomic DNA <1 µg, plasmid or viral DNA 1-50 ng, cDNA synthesis reaction 1-5 µL (or <100 ng).

3. Dispense the master mix into wells of PCR plate or tubes.
4. Add in each well the DNA sample. Mix well by pipetting and spin down the samples.
5. Seal the PCR plate with optical film or tube caps.
6. Set up thermocycler with your optimized program conditions (if applicable, select fast mode on the instrument*):

Suggested thermal cycling conditions:

Process	Cycles	Temperature	Time
Enzyme Activation	1 x	95 °C	2-3 min
Denaturation	40 x	95 °C	5 sec
Annealing/ Extension		60-65°C	20-30sec**
Dissociation/Melt Analysis	See your instrument guidelines for setup.		

*For efficient qPCR, under fast cycling conditions, it is recommended to amplify DNA fragments ranging from 80-200bp.

** Select the shortest time possible but not less than 20 sec and do not exceed 30 seconds. Primers should have melting temperature of approximately 60°C.

In case of 3-step cycling, anneal at optimal annealing temperature for 20 sec and minimum time necessary at 72°C for data acquisition (according to manufacturer's guidelines).

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