Tiaris™ Tissue Direct PCR

Ordering Info

TBK1015, 80 reactions of 25 µL

TBK1016, 480 reactions of 25 µL

Description

Tiaris™ Tissue Direct PCR kit is designed to perform PCR directly from tissue samples with no prior DNA purification. The kit combines a simple but efficient DNA extraction method with direct amplification using our Fast Hot Start DNA polymerase in a convenient and easy-to-use manner. This kit can be used with a variety of samples, including mouse ear and tail, FFPE tissue and human hair. DNA extraction is carried out in a single tube, without the need of multiple washing steps, therefore minimizing the risk of contamination. The Fast Hot Start DNA polymerase Master Mix contains premixed gel loading dye which allows direct sample loading on agarose gel.

Features

- DNA is ready for PCR in only 15 minutes.
- Amplicon size: up to 5kb.
- Extension Rate 4-8 kb/min.
- Compatible with fast and standard PCR program.

Kit Components

Components	TBK1015	TBK1016
Lysis Buffer A	1.6 mL	6 x 1.5 mL
Lysis Buffer B	800 µL	6 x 800 µL
PCR-grade water	1 mL	6 x 1 mL
Fast Hot-Start DNA polymerase Master Mix(2x)	1 mL	6 x 1 mL


Order Info Kit Components: Lysis Buffer A (TBK1015-1) | Lysis Buffer B (TBK1015-2) | Fast Hot-Start DNA polymerase Master Mix, 2x (TBK1028) | PCR Grade Water (TBB0303).

Storage

Shipped on blue ice. Upon receipt, kit components should be immediately stored at -20°C. No loss of performance is detected after 20 freeze/thaw cycles.

Applications

- Direct PCR
- Mouse Genotyping
- Genetic Screening
- Multiplex PCR

PROTOCOL

A. EXTRACTION

1. Prepare for each sample **100 µL reaction mixture** consisting of **70 µL PCR-grade water**, **20 µL Lysis Buffer A**, and **10 µL Lysis Buffer B**.
2. Cut 3-30 mg tissue, ear punch, or mouse tail sample and place into a 1.5 mL microcentrifuge tube.
3. Cover each tissue sample (see Part C) with 100 µL of **reaction mixture**. Mix thoroughly by pipetting or vortexing.
4. Incubate at 75°C for 5 minutes (*lysis and protein denaturation*).
5. Heat-inactivate by incubation at 95°C for 10 minutes.
6. Centrifuge at 14.000-16.000g for 1 minute to pellet debris and transfer the cleared supernatant to a new DNase-free microtube. Any remaining visible residue will not have a negative impact on the downstream PCR.
7. Store at -20°C for up to several months, or proceed with PCR directly.

B. PCR AMPLIFICATION

1. Gently vortex and briefly centrifuge Fast Hot Start DNA polymerase Master Mix after thawing.
2. Place a tube on ice and add the following components for each 25 µL reaction. Prepare sufficient master mix for the number of reactions. Consider one or two extras:

Components	Volume	Final Concentration
Fast Hot Start DNA polymerase Master Mix	12.5 µL	1x
Forward primer (5 pmol/µL)	2 µL	400 nM
Reverse primer (5 pmol/µL)	2 µL	400 nM
Diluted DNA extract (Step 5)	1 µL	*
PCR Grade Water, nuclease-free	7.5 µL	
Final Volume	25 µL	

* Please note that it might be necessary to optimize the amount of template. We recommend to dilute the DNA extract 10-fold with PCR-grade water (e.g. 5 µL of DNA extract + 45 µL PCR Grade Water) and then use, as indicated in the table, 1 µL as template. This 10-fold dilution ensures the dilution of PCR inhibitors. If the DNA concentration is too low, one can dilute less (e.g. 5-fold or undiluted), however, this might result in a less efficient amplification as the concentration of inhibitors will also be higher and it might be necessary to adapt the cycling conditions (number of cycles and/or extension time). Alternatively, results might be improved if template DNA is diluted even more (e.g. 20-fold).

3. Dispense the master mix into wells of PCR plate.
4. Gently vortex and spin down the samples.
5. Add in each well the DNA sample (*Diluted DNA extract**). Mix well by pipetting.
6. Seal the PCR plate with optical film.

7. Set-up PCR amplification. The amplification parameters should be optimized for individual primers, template, and thermal cycler.

Step	Cycles	Temperature	Time
Initial denaturation	1 x	95 °C	3:00
Denaturation		95 °C	0:15
Annealing	40 x	T _m	0:15
Extension		70 °C	**
Final Extension	1 x	70 °C	3:00
Conservation	1 x	4 °C	∞

** 2 to 75 seconds at 70°C for extension (2 seconds for targets below 1 kb and 15 seconds per kb for target DNA up to 5kb. 90 seconds total per cycle in case of multiplex PCR, independent of the sizes of the amplicons).

C. SAMPLE GUIDELINES

Sample (fresh or frozen)	Amount	Extraction Volume
Mouse tail	1 - 2 mm (~ 5 mg)	100 µL
Mouse ear	2 - 4 mm ² (~ 5 mg)	100 µL
Mammalian tissue	3 - 30 mg	100 µL
Hair follicle	1 - 10 follicles	100 µL

Sample amounts can be slightly increased for better yields, but too much material may cause inefficient lysis and PCR inhibition.

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