Tiaris™ Blood Direct PCR Kit

Ordering Info

TBK1018, 80 reactions of 25 µL

TBK1019, 480 reactions of 25 µL

Description

Tiaris™ Blood Direct PCR kit allows direct PCR analysis of whole blood sample, without the need of DNA purification and sample preparation. The kit combines a simple but efficient DNA extraction method with direct amplification using our Fast Hot Start DNA polymerase in a convenient and easy-to-use manner. This kit is applicable for fresh blood, blood stored at 4°C and dried blood spot stored on DBS commercial cards. DNA extraction is carried out in a single tube, without the need of multiple washing steps, therefore minimizing the risk of contamination. The Fast Hot Start DNA polymerase Master Mix contains premixed gel loading dye which allows direct sample loading on agarose gel.

Features

- DNA is ready for PCR in only 15 minutes.
- Amplicon size: up to 5kb.
- Extension Rate 4-8 kb/min.
- Compatible with fast and standard PCR program.

Kit Components

Components	TBK1018	TBK1019
Lysis Buffer A	1.6 mL	6 x 1.6 mL
Lysis Buffer B	800 µL	6 x 800 µL
PCR-grade water	1 mL	6 x 1 mL
Fast Hot Start DNA polymerase Master Mix(2x)	1 mL	6 x 1 mL


Order Info Kit Components: Lysis Buffer A (TBK1018-1) | Lysis Buffer B (TBK1018-2) | Fast Hot Start DNA polymerase Master Mix, 2x (TBK1028) | PCR Grade Water (TBB0303).

Storage

Shipped on blue ice. Upon receipt, kit components should be immediately stored at -20°C. No loss of performance is detected after 20 freeze/thaw cycles.

Applications

- Direct PCR
- Genotyping (such as gene detection) analysis
- Multiplex PCR

PROTOCOL

A. EXTRACTION

1. Prepare for each sample **100 µL** reaction mixture consisting of **70 µL PCR-grade water**, **20 µL Lysis Buffer A**, and **10 µL Lysis Buffer B**.
2. Cover each sample (*see Part C*) with 100 µL of reaction mixture and incubate at 75°C for 5 minutes (*lysis and protein denaturation*).
3. Heat-inactivate by incubation at 95 °C for 10 minutes.
4. Centrifuge at 14,000-16,000g for 1 minute to pellet debris and transfer the cleared supernatant to a new DNase-free microtube. Any remaining visible residue will not have a negative impact on the downstream PCR.
5. Store at -20°C for up to several months, or proceed with PCR directly.

B. PCR AMPLIFICATION

1. Gently vortex and briefly centrifuge Fast Hot Start DNA polymerase Master Mix after thawing.
2. Place a tube on ice and add the following components for each 25 µL reaction. Prepare sufficient master mix for the number of reactions. Consider one or two extras:

Components	Volume	Final Concentration
Fast Hot Start DNA polymerase Master Mix	12.5 µL	1X
Forward primer (5 pmol/µL)	2 µL	400 nM
Reverse primer (5 pmol/µL)	2 µL	400 nM
Diluted DNA extract (<i>Step 5</i>)	1 µL	*
PCR Grade Water, nuclease-free	7.5 µL	
Final Volume	25 µL	

* Please note that it might be necessary to optimize the amount of template. We recommend to dilute the DNA extract 10-fold with PCR-grade water (e.g. 5 µL of DNA extract + 45 µL PCR Grade Water) and then use, as indicated in the table, 1 µL as template. This 10-fold dilution ensures the dilution of PCR inhibitors. If the DNA concentration is too low, one can dilute less (e.g. 5-fold or undiluted), however, this might result in a less efficient amplification as the concentration of inhibitors will also be higher and it might be necessary to adapt the cycling conditions (number of cycles and/or extension time). Alternatively, results might be improved if template DNA is diluted even more (e.g. 20-fold).

3. Dispense the reaction mix into wells of PCR plate.
4. Gently vortex and spin down the samples.
5. Add in each well the DNA sample (*diluted DNA extract**). Mix well by pipetting.
6. Seal the PCR plate with optical film.

7. Set-up PCR amplification. The amplification parameters should be optimized for individual primers, template, and thermal cycler.

Step	Cycles	Temperature	Time
Initial denaturation	1 x	95 °C	3:00
Denaturation		95 °C	0:15
Annealing	40 x	T _m	0:15
Extension		70 °C	**
Final Extension	1 x	70 °C	3:00
Conservation	1 x	4 °C	∞

** 2 to 75 seconds at 70°C for extension (2 seconds for targets below 1 kb and 15 seconds per kb for target DNA up to 5kb. 90 seconds total per cycle in case of multiplex PCR, independent of the sizes of the amplicons).

C. SAMPLE GUIDELINES

Sample (fresh or frozen)	Amount	Extraction Volume
Blood (mammalian)	2 – 8 µL (with EDTA as anti-coagulant)	100 µL
Dried Blood Spot Cards	2 mm ²	100 µL

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