IPTG

(Isopropyl- β -D-thiogalactopyranoside)

Ordering info

TBR0114, IPTG 5g

TBR0115, IPTG Solution 1M, 5x 1.5 mL

TBR0116, IPTG Solution 0.1M, 10x 1.5 mL

Description

IPTG is a synthetic compound widely used in protein expression driven by lac promoter. This compound is a non-metabolizing analogous of lactose that triggers transcription of the lac promoter.

CAS: 367-93-1

Formula: $C_9H_{18}O_5S$

MW: 238,30 g/mol

Storage

Store at $-20^{\circ}C$. Stable for 3 years.

Quality Control

- Purity is tested by HPLC (>99%).

Features

- Highest purity (>99%).
- White to yellow powder with a sulfur odor.
- Dioxane free.
- Solubility ($20^{\circ}C$): > 20% in water.

References

Sambrook, J. and Russell, D.W. (2001) Molecular cloning: a laboratory manual, 3rd ed., Cold Spring Harbor Laboratory Press.

Applications

- Induction of protein expression.
- In combination with X-GAL is used to prepare X-GAL blue/ white selective plates.

Also available:

X-GAL (TBR0111, 1 g)

Ampicillin (TBR0112, TBR0113)

PROTOCOL

X-GAL plates includes X-GAL, IPTG and an antibiotic resuspended in a culture media, generally Luria broth agar. X-GAL plates are usually prepared in two different ways.

Using LB Agar media melted

1. Melt LB Agar media.
2. Cool until 50°C and add the following components (volume indicated for 100 mL LB Agar)

Components	Volume	Final Concentration
X-GAL Solution 2% (TBR0111)	250 µL	0.5 mg/mL
IPTG 0,1 M	134 µL	134 µM
Antibiotic	*	*
LB Agar Media	100 mL	1x

* Depends of each antibiotic.

3. Mix well and pour the mixture on sterile petri dishes.
4. Protect X-GAL plates from light.

Using LB Agar + Antibiotic plates

1. On the surface of the plate, add 40 µL X-GAL Solution 2% and 40 µL IPTG 0.1 M.
2. Spread on the plate with a sterile Drigalsky spatula. Allow the liquid absorption (~ 20 minutes) prior bacterial plating.

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