# TIARISTAIN™ Green Safe

## (Molecular Biology Grade)

### Ordering info

TBR0225, TIARISTAIN™ Green Safe, 20,000x, 20 µL

TBR0226, TIARISTAIN™ Green Safe, 20,000x, 1 mL

### Description

TIARISTAIN™ Green Safe is a new formulation to stain nucleic acids in agarose and polyacrylamide gels. It emits a green fluorescence when bound to nucleic acids. TIARISTAIN™ Green Safe has a strong excitation peak at 485 nm enabling the use of less DNA destructive blue light instead UV light (better low pressure lamp, 254 nm). The dye has also two excitation peaks in the ultraviolet range, at 270 nm and 290 nm). The maximum emission is at 525 nm.

### Features

- Non carcinogenic high quality product.
- Excited with both blue or ultraviolet light and has a maximum emission at 525 nm (green).
- High sensitivity, 0.2-0.6 ng DNA per band.
- Flexible use, valid for pre-casting and post-staining use.

### Applications

- dsDNA, ssDNA and RNA staining in agarose gels.
- Non-carcinogenic alternative to ethidium bromide nucleic acid stain.

### Storage

Store at 2-8 °C. Protect from light.

The product is shipped at room temperature.

### Quality Control

DNase/ RNase activity not detected.

### Recommended devices

Transilluminators with blue light (470-520 nm), ultraviolet light (preferable with low pressure lamp 254 nm; high pressure 365 nm but excitation is less efficient).

Also available:

**INTENSE™ INDIGO™ Loading Buffer 6x (TBB0325)**

**COBALT™ Loading Buffer 6x (TBB0321)**

**GREENY™ Loading Buffer 6x (TBB0322)**

**CORAL™ Loading Buffer 6x (TBB0323)**

**COLLAGE™ Loading Buffer 6x (TBB0324)**

## PROTOCOL

; Wear gloves while handling!

### I. PRE-CASTING

1. Prepare an agarose gel (0.5 - 3% in TAE or TBE Buffer 1x).
2. Heat the solution until all the agarose is completely melted.
3. Add TIARISTAIN™ Green Safe and mix thoroughly,

V (Agarose Gel)	V (TIARISTAIN™ Green Safe, 20,000x)
100 mL	5 µL
50 mL	3 µL

4. Cool the solution and pour into the tray.
5. Once the gel is solid, load the samples.
6. For best results, add 1 µL /10 mL running buffer near the anode.

### II. POST STAINING

*Post-staining gives a better sensitivity than precast gels.*

1. Prepare the Post-Staining Solution adding **20 µL TIARISTAIN™ Green Safe 20,000x** to 100 mL TAE or TBE Buffer 1x (agarose gel thickness <0.5 cm). Mix well.

*Post-Staining Solution could be re-used 2-3 times. Store in the dark at room temperature.*

2. Soak the gel at room temperature for 10-60 minutes.

*Incubation time depends on agarose percentage and gel thickness.*

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