L-Alanyl-L-Glutamine 100x

(Cell Culture Grade)

Ordering info

TBR0322, L-Alanyl-L-Glutamine 100x, 100 mL

Description

L-Alanyl-L-Glutamine 100x is a concentrated sterile solution of **stable L-glutamine**, used as a supplement in cell culture medium. The dipeptide is metabolized within the cells to yield L-Glutamine plus the second amino acid. This results in a more consistent delivery of L-Glutamine to cells and avoids toxic build-up of ammonia in cell cultures. This feature can be especially important for ammonia-sensitive cell lines.

L-Alanyl-L-Glutamine prevents the intramolecular cyclization reaction associated with solutions of L-Glutamine, allowing the formulation of cell culture media containing L-Glutamine that may be stored at 4°C for extended periods. Solutions containing these derivatives can even be autoclaved without appreciable degradation of the product (30 minutes at 121°C results in <5% loss of the product).

L-glutamine is an essential amino acid that plays a key role in cellular metabolism, serving as a major source of nitrogen and carbon, which are essential for synthesizing proteins, nucleotides, and other molecules required for cell proliferation.

The concentration of L-glutamine needed varies depending on the cell type and medium formulation. Classical cell culture media like DMEM, GMEM, and IMDM typically contain L-glutamine concentrations ranging from 2 mM to 4 mM, which are adequate for most cell lines. However, certain cells, such as hybridoma cells, may require higher concentrations for optimal growth. In specialized media, such as serum-free or protein-free formulations, L-glutamine is often supplemented at levels around 2.5 mM to 6 mM.

Features

- Sterile colorless solution.
- pH = 7.0 ± 0.5
- Endotoxin free (< 10 EU/ mL)
- Composition (100x): L-Alanyl-L-Glutamine 200 mM.

Storage

Store at -20 °C.

For an appropriate use of the product, prepare aliquots and freeze at -20°C. Avoid repeated freeze/thaw cycles.

Applications

Support the growth and maintenance of cells.

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