Trypsin-EDTA

In DPBS

(Cell Culture Grade)

Ordering info

TBZ0342, Trypsin-EDTA 1x, 100 mL

TBZ0344, Trypsin-EDTA 10x, 100 mL

Description

Trypsin-EDTA is a well-known solution used as standard cell detachment solution to detach cells from standard tissue culture plastic ware and adhesion coated plastic ware. It contains pancreatic porcine trypsin at 0.05% (1x) or 0.5% (10x) in DPBS pH=7.3 ± 0.3. Trypsin is a serine protease that hydrolyzes proteins at the carboxyl side of the Lysine or Arginine. EDTA presence enhances trypsin action, and reduces the required trypsin concentration for effective hydrolysis. The solution does not include calcium and magnesium.

Storage

Store at -20°C. Stable 2 years. Defrost at 4°C !

The product is shipped on blue ice.

Material required (not included)

- PBS 1x pH=7.4 (TBB0360, TBB0361)
- PBS 1x pH=7.2 (TBB0366, TBB0367)
- Specific Growth Media
- Serum
- Trypan Blue Cell Viability Indicator (TBB0402, TBB0403)

Features

- Sterile (by filtration).
- Ready to use.
- Composition 10x: 0.5% Porcine Trypsin in Dulbecco' PBS pH=7.3 ± 0.3.
- Ca²⁺ and Mg²⁺ free solution.
- Colorless solution.

Applications

- Dissociating adherent cells from surfaces.
- Routine cell passage.
- Create single cell suspension for accurate cell counting.

PROTOCOL

1. Prewarm **Trypsin-EDTA 1x** and **PBS 1x** at 37 °C.
2. Carefully, aspirate all the culture media off.
3. Wash the cells with **PBS 1x**.

The cells should be thoroughly covered with PBS 1x. This rinse is instantaneous but the solution can remain on the cells for up to 4 hours.

4. Pour off the PBS rinsing solution.
5. Add **Trypsin-EDTA 1x** to each vessel as follows:

Flasks	Volume
25 cm ²	0.5 - 1 mL
75 cm ²	2.5 - 5 mL
150 cm ²	5- 10 mL

6. Incubate at 37 °C for 2-3 minutes*.

** Check cell morphology under microscope to verify if cell have detached; they look rounded and refractile. If not, continue incubation for another 2 minutes and check the cells again. Avoid exposures longer than 10 minutes. The time depends on cell type, density, culture age, etc. Ever exceed 10 minutes.*

7. Once cells are detached, hit the flask against the palm of the hand to ensure all cells have been dislodged.
8. To inhibit trypsin activity, add prewarmed **Growth media containing serum**. Add at least 10 times the volume of Trypsin-EDTA used.
For serum free cultures, inhibit trypsin activity adding an equimolar concentration of soybean trypsin inhibitor.
9. By gentle pipetting, recover the cells in a fresh tube.
10. Centrifuge at 100-300 g for 5-10 minutes and discard the supernatant.
11. Resuspend the cell pellet with prewarmed growth medium.
12. Determine cell viability using **Trypan Blue Cell Viability Indicator** (TBB0402-0403).

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