Blood Genomic DNA Purification Midi Kit

(for processing 2 mL blood sample)

Ordering info

TBK0129, 20 reactions

TBK0130, 200 reactions

Description

Blood Genomic DNA Purification Midi Kit is an optimized kit to obtain high molecular weight genomic DNA from blood. It is based on salting-out principle to produce higher quantity and quality of DNA. Suitable for EDTA, sodium citrate and heparin-anticoagulated fresh or frozen blood.

Features

- High yield and purity, 3-6 µg, A₂₆₀/A₂₈₀ ~1.8.
- Scalable, easily to process many samples simultaneously.
- No phenol extraction.
- Fast and easy protocol.
- Cost-effective.

Applications

DNA obtained is suitable for downstream molecular biology applications such as PCR, enzymatic digestion for cloning or Southern, genotyping, etc.

Quality Control

DNA isolation from 200 µL whole human blood is checked by: integrity (agarose gel electrophoresis), quantity and quality (A₂₆₀/A₂₈₀ = 1.8 ± 0.2).

Kit Components

Components	TBK0129	TBK0130
BG1 Buffer	125 mL	1250 mL
BG2 Buffer	50 mL	430 mL
BG3 Buffer	20 mL	150 mL
Proteinase K*	2x30 mg ^a	500 mg ^b
Elution Buffer	15 mL	45 mL

Order Info Kit Components: BG1 Buffer (TBB0501) | BG2 Buffer (TBB0502) | BG3 Buffer (TBB0503) | Proteinase K (TBZ0305, TBZ0308) | Elution Buffer (TBB0510).

* To prepare a 20 mg/mL solution add:

^a 1.5 mL Water (Molecular Biology Grade) to 30 mg Proteinase K powder.

^b 25 mL Water (Molecular Biology Grade) to 500 mg Proteinase K powder.

Store Proteinase K solution in aliquots at -20°C.

Storage

Store the kit at 25°C.

Store Proteinase K at 2-8°C (short storage) or -20°C (long storage).

Material required (not supplied)


- Isopropyl Alcohol (CAS 67-63-0).
- Ethanol 70%.
- 15 mL tubes.

PROTOCOL

1. Transfer **2 mL of blood** to a 15 mL tube.
2. Add **6 mL BG1 Buffer** and mix by vortex
3. Incubate 10 minutes at room temperature.
4. Centrifuge at 3,000 g for 2 minutes and remove the supernatant leaving a residual volume of ~200 μ L.
5. Vortex 10-15 seconds until pellet are resuspended.
6. Add **100 μ L Proteinase K** and **2 mL BG2 Buffer**.
7. Pipet up and down gently to mix.
8. Incubate at 55°C for 0.5-1 hour.
9. Add **700 μ L BG3 Buffer** and mix vigorously by vortex.
10. Centrifuge at 3,000 g for 5 minutes.

If no dark brown pellet is seen, incubate the sample on ice 5 minutes and repeat step 9.

11. Recover the supernatant by transferring it to a fresh 15 mL tube.
12. Add **2 mL Isopropanol**.
13. Mix by gentle inversion (~ 50 times).
14. Centrifuge at 3,000 g for 3-5 minutes and discard the supernatant with a pipette.
15. Add **2 mL Ethanol 70%**.
16. Centrifuge at 3,000 g for 2 minutes and remove the supernatant using a pipette.
17. Invert and dry the tubes on clean paper towels for 5-10 minutes.
18. Add **200 μ L Elution Buffer** to resuspend the DNA.
19. Check DNA integrity on agarose electrophoresis gel and quantity and quality by spectrophotometry.
20. Store at -20°C.

Blood Genomic DNA Purification Midi Kit

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