

## M-MuLV Reverse Transcriptase

### Ordering Info

TBZ0321, 20,000 U (200 U/  $\mu$ L)

TBZ0322, 100,000 U (200 U/  $\mu$ L)

### Description

M-MuLV Reverse Transcriptase is a recombinant reverse transcriptase from Moloney Murine Leukemia Virus produced in *Escherichia coli*. Highly purified the enzyme is an useful RNA-dependent DNA polymerase to synthesize cDNA using a RNA template and an oligo(dT) primer or a specific reverse primer.

CAS: 9068-38-6

MW: 69 kDa

### Features

- The enzyme lacks of 3'→5' and 5'→3' exonuclease activity.
- The enzyme has not RNase H activity.

### Applications

- First-strand cDNA synthesis experiments
- RT-PCR
- RT-qPCR

### Quality Control

Functionally tested in a quantitative RT-PCR of actin gene using as template dilutions of total RNA transcribed using 200 units of enzyme.

### Kit Components

Components	TBZ0321	TBZ0322
M-MuLV Reverse Transcriptase (200 U/ $\mu$ L)	100 $\mu$ L	500 $\mu$ L
M-MuLV Reaction Buffer (10x)	1.5 mL	2 x 1.5 mL

**Order Info Kit Components:** M-MuLV Reaction Buffer (10x) (TBB0317).

### Storage

Store at -20°C. Shipped in blue ice.

### Unit Definition

One unit is defined as the amount of enzyme required to incorporates 1 nmol of dTTP into acid-insoluble material in 10 min at 37 °C using poly(A) : oligo (dT) as template : primer.

### Material required (not supplied)

- dNTPs mix 2.5 mM each (TBR0209)
- PCR Grade Water (TBB0303)
- Ribonuclease Inhibitor
- PCR Tubes
- Oligo(dT) and specific primers

## PROTOCOL

1. Thawing all components on ice. Vortex them and centrifugate.
2. On ice, prepare a mix of the following components,

Reaction Components	Final Concentration	Volume
RNA template	Total RNA: 10 pg - 5 µg mRNA: 10 pg - 0.5 µg	x µL
dNTPs Mix 40 mM TOTAL	0.5 mM	1 µL
Oligo(dT) <sub>12-18</sub> (50-60 µM) or Specific Reverse Primer (2 µM) or Random Hexamer (50-250 ng/ µL)	2.5 – 3 µM 0.1 µM 2.5 – 12.5 ng/ µL	1 µL
Nuclease free Water		up 16 µL
<b>Final Volume</b>		<b>16 µL</b>

3. [Optional] Incubate at 65 °C, 5 minutes.  
*This step is recommended for sequences with high GC content.*
4. Spin the tubes and add the following components on ice:

Reaction Components	Final Concentration	Volume
M-MuLV Reaction Buffer (10x)	1 x	2 µL
Ribonuclease Inhibitor		1 µL
M-MuLV Reverse Transcriptase (200 U/ µL)	10 U/ µL	1 µL
<b>Final Volume</b>		<b>20 µL</b>

5. Incubate at 37 °C, 50 minutes.  
*If you are using random hexamers, firstly incubate at 25°C, 15 minutes and then at 37 °C, 50 minutes.*
6. To inactivate the enzyme, incubate at 70°C, 15 minutes.
7. Store the synthesized cDNA at -20°C.