STOUT™-Green PCR Master Mix (2x)

Ordering info

TBK0031, 5 reactions (sample)

TBK0032, 4 x 1.25 mL (200 reactions)

TBK0033, 8 x 1.25 mL (400 reactions)

Description

STOUT™-Green PCR Master Mix (2x) is an optimized formulation to enable successful amplification in a process with reduced manipulation. The formula contains MgCl₂, pure dNTPs, PCR enhancer and our highly purified thermostable STOUT™ Recombinant Taq DNA Polymerase. The master mix is a dense solution for direct loading of PCR products onto agarose gel and it includes two tracking dyes to monitor the sample electrophoresis.

Features

- **Ready to use**, avoiding mistakes in PCR reaction preparation.
- **Optima formulation** in routine PCR, requiring to add only specific primers and template.
- **High Thermostability**, the enzyme half-life at 94 °C is 40 minutes.
- **TA Cloning**, addition of 3'adenine at the end of PCR fragment.
- Error Rate 1–20x10⁻⁵ errors/bp per cycle.

Applications

- Daily research: standard PCR for targets up 5kb, recombinant clone or cell line checking, etc.
- Generation of PCR fragments for TA cloning.

Kit Components

Components		TBK0032	TBK0033
STOUT™-Green Master Mix (2x)	PCR	4 x 1.25 mL	8 x 1.25 mL
PCR Grade Water, nuclease free		3 x 2 mL	5 x 2 mL

Order Info Kit Components: PCR Grade Water, nuclease free (TBB0303).

Storage

Store at -20°C. Shipped in blue ice.

Unit Definition

One unit of STOUT™ Recombinant Taq DNA polymerase catalyzes the incorporation of 10 nanomoles of dNTPs into acid-insoluble material in 30 minutes at 74 °C.

Quality Control

Functionally tested in a 1kb PCR amplification (GC 52%).

Material required (not supplied)

- PCR Tubes.
- Primers.

PROTOCOL

I. PREPARING REACTIONS

1. Gently vortex and briefly centrifuge STOUT™-Green PCR Master Mix (2x).
2. Place a tube on ice and add the following components for each 50 µL reaction. Prepare sufficient master mix for the number of reactions. Consider one or two extras:

Components	Volume	Final Concentration
STOUT™-Green PCR Master Mix (2x)	25 µL	1x
Forward primer 15 µM (15 pmol/µL)	1 µL	0.75 µM
Reverse primer 15 µM (15 pmol/µL)	1 µL	0.75 µM
PCR Grade Water, nuclease-free	up 50 µL	
Template DNA		**
Final Volume	50 µL	

** Optimal amounts of template DNA are 10-30 ng for both plasmid and phage DNA, and 0.3 -1 µg for genomic DNA.

3. Gently vortex and spin down the mix. Aliquot the mix into individual PCR tubes.
4. Add template DNA.

II. PCR SETUP

1. Perform PCR using recommended thermal cycling conditions:

Step	Cycles	Temperature	Time
Initial denaturation	1 x	94 °C	5:00
Denaturation		94 °C	0:35
Annealing	25 - 30 x	Tm	0:35
Extension		72 °C	1:00 per kb
Final Extension	1 x	72 °C	10:00
Conservation	1 x	4 °C	∞

2. Store the samples at -20°C until use.
3. Load the samples on an agarose gel.

In an 0.7% agarose gel, blue dye migrates at 8,000 bp while yellow dye migrates at 100 bp in TAE 1x buffer.

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