

High-Q™-Spin-Column Tissue Genomic DNA Purification Kit

Ordering info

TBK0162. 3 reactions (sample)

TBK0163. 50 reactions

TBK0164. 200 reactions

Description

High-Q™-Spin-Column Tissue Genomic DNA Purification Kit is an easy silica-membrane-based system for genomic DNA purification from a wide variety of tissue samples. An optimized lysis buffer guarantees a good yield while the use of High-Q™ Spin Columns allow a good quality DNA, suitable for downstream applications.

Features

- Starting material up to 100 mg of tissue sample.
- Typical yields are **0.5- 50 µg of DNA** depending on the material tissue used.
- No organic extraction, no ethanol precipitation.
- High DNA purity, the isolated DNA is ready to use for downstream molecular biology applications.
- Easy and Fast protocol.

Applications

- Purification of DNA from animal and human tissue.
- Purification of DNA tissue using different starting materials: frozen, fresh or stabilized with Tiaris™ Tissue Protect (TBB0431).
- DNA obtained is suitable for downstream molecular biology applications such as PCR, enzymatic digestion for cloning or Southern, genotyping, etc.

Quality Control

DNA purified is checked by: integrity (agarose gel electrophoresis), quantity and quality ($A_{260}/_{280} = 1.9 \pm 0.2$).

Kit Components

Components	TBK0163	TBK0164
High-Q™ Spin Column with Collection Tubes	50	200
Tissue-1 Buffer	30 mL	110 mL
Tissue-2 Buffer	12 mL ^a	40 mL ^b
WB2 Buffer	12 mL ^c	45 mL ^d
Elution Buffer	15 mL	25 mL
Proteinase K	2x 30 mg ^e	6x 30 mg ^e
Proteinase K Resuspension Buffer	2x 1.5 mL	6x 1.5 mL
RNase A	15 mg ^f	3 x 15 mg ^f
RNase A Resuspension Buffer	1 x 1.5 mL	3 x 1.5 mL

Order Info Kit Components: High-Q™ Spin Column with Collection Tubes (TBM0010) | Tissue-1 Buffer (TBB0579) | Tissue-2 Buffer (TBB0580) | WB2 Buffer (TBB0512) | Elution Buffer (TBB0510) | Proteinase K (TBZ0303) | Proteinase K Resuspension Buffer (TBB0546) | RNase A (TBZ0318) | RNase A Resuspension Buffer (TBB0309).

Components for samples are ready to use!

Before its use:

- ^a Add 18 mL isopropanol and mix well.
- ^b Add 60 mL isopropanol and mix well.
- ^c Add 48 mL absolute ethanol and mix well
- ^d Add 180 mL absolute ethanol and mix well.
- ^e Add 1.5 mL Proteinase K Resuspension Buffer and mix well. Store at -20°C.
- ^f Add 1.5 mL RNase A Resuspension Buffer and mix well. Store at -20°C.

Storage

Store the kit at 25°C.

Store Proteinase K and RNase solutions at -20°C.

Material required (not supplied)

- 1.5 mL Tubes.
- Ethanol (CAS 64-17-5).
- Isopropanol (CAS 67-63-0).

PROTOCOL

- Grind up to 50-100 mg of tissue material in liquid nitrogen using a mortar and a pestle. With a freeze spatula, collect the powder into a 1.5 mL tube. Commercially available equipment for homogenization also can be used.
- Add **500 µL Tissue-1 Buffer** and mix by vortex briefly.
- Add **40 µL Proteinase K (20 mg/mL)**. Mix well.
- Incubate at 65°C until the tissue is completely lysed. Lysis is usually complete in 1-3 h. Some samples as rodent tail and ears should be lysed overnight without affect the procedure results. Mix 2-3 times by inversion during incubation.
- Centrifuge at 13,000 g for 5 minutes. Transfer the supernatant to a fresh tube.
- [Optional]** Add **20 µL RNase-A 20 µg/ µL** and incubate for 5 minutes at room temperature.
- Add **400 µL Tissue-2 Buffer** and mix by inversion.
Check isopropanol has been added to Tissue-2 Buffer (✓).
- Transfer up **700 µL mixture** to a High-Q™ Spin Column placed into a Collection Tube.
- Centrifuge at 10,000 g for 1 minute. Remove the flow-through and place back the High-Q™ Spin Column into a Collection Tube. Repeat steps 8 and 9 with the remaining mixture.
- Add **500 µL WB2 Buffer**.
Check absolute ethanol has been added to WB2 Buffer (✓).
- Centrifuge at 10,000 g for 30 seconds. Discard the flow-through and place back the High-Q™ Spin Column into a Collection Tube. Repeat steps 10 and 11 one more time.
- To dry High-Q™ Spin Column and eliminate residual ethanol, centrifuge again at 13,000 g for 1 minute.
- Place the High-Q™ Spin Column into a clean 1.5 mL Tube. Add **50-100 µL prewarmed Elution Buffer or Water (Molecular Biology Grade)**.
Prewarm Elution Buffer or Water at 70°C.
- Incubate at room temperature, 2 minutes. Centrifuge at 13,000 g for 1 minute to elute purified DNA.
- Check DNA quantity by spectrophotometry and quality on agarose electrophoresis gel.
- Store DNA at -20°C.

