

High-Q™-Spin-Column Fungal Genomic DNA Purification Kit

Ordering info

TBK0254, 3 reactions (sample)

TBK0255, 50 reactions

TBK0256, 200 reactions

TBK0257, 250 reactions

Description

High-Q™-Spin-Column Fungal Genomic DNA Purification Kit is an easy silica-membrane-based system for DNA purification from a wide variety of fungal species. An optimized lysis buffer guarantees a good yield while the use of High-Q™ Spin Columns allow a good quality DNA, suitable for downstream applications.

Features

- Easy and Fast protocol.
- Typical yields are 2- 50 µg of genomic DNA.
- No organic extraction, DNA is isolated without phenol, chloroform or ethanol precipitation.
- High DNA purity, the isolated DNA is ready to use for sensitive downstream molecular biology applications.
- Isolation of DNA from a wide range of samples.

Applications

- Purification of DNA from fungal samples.
- DNA obtained is suitable for downstream molecular biology applications such as PCR, enzymatic digestion for cloning or Southern, genotyping, etc.

Quality Control

DNA purified is checked by: integrity (agarose gel electrophoresis), quantity and quality ($A_{260}/_{280} = 1.8 \pm 0.2$).

Kit Components

Components	TBK0255	TBK0256	TBK0257
High-Q™ Spin Column with Collection Tubes	50	200	250
BFL1 Buffer	25 mL	85 mL	110 mL
BFL1 Additive	3 mL	10 mL	12 mL
BFL2 Buffer	1.5 mL ^a	5 mL ^b	6 mL ^c
BFL3 Buffer	7 mL ^d	28 mL ^e	35 mL ^f
BFL4 Buffer	12 mL ^g	45 mL ^h	2 x 30 mL ⁱ
Elution Buffer	15 mL	25 mL	25 mL
Proteinase K	1.1 mL [*]	3 x 30 mg ^{**}	4 x 30 mg ^{**}
Proteinase K Buffer	-	3 x 1.5 mL	4 x 1.5 mL

Order Info Kit Components: High-Q™ Spin Columns (TBM0010) | BFL1 Buffer (TBB0547) | BFL1 Additive (TBB0548) | BFL2 Buffer (TBB0549) | BFL3 Buffer (TBB0550) | BFL4 Buffer (TBB0551) | Elution Buffer (TBB0510) | Proteinase K 20 mg/ mL (TBZ0308)* | Proteinase K 30 mg (TBZ0305)** | Proteinase K Buffer (TBB0546).

Components for samples are ready to use!

Before its use:

- ^a Add 13.5 mL isopropanol and mix well.
- ^b Add 45 mL isopropanol and mix well.
- ^c Add 54 mL isopropanol and mix well.
- ^d Add 21 mL absolute ethanol and mix well.
- ^e Add 84 mL absolute ethanol and mix well.
- ^f Add 105 mL absolute ethanol and mix well.
- ^g Add 48 mL absolute ethanol and mix well.
- ^h Add 180 mL absolute ethanol and mix well.
- ⁱ Add 120 mL absolute ethanol to each bottle. Mix well.

Storage

Store the kit at 25°C and Proteinase K at -20°C.

Material required (not supplied)

- 1.5 mL Tubes.
- Ethanol (CAS 64-17-5).
- Isopropanol (CAS 67-63-0).

PROTOCOL

- Grind up to 100 mg of fungal mycelia in liquid nitrogen using a mortar and a pestle. With a freeze spatula, collect the powder into a 1.5mL tube. Commercially available equipment for homogenization can also be used.
- Add **400 µL BFL1 Buffer** and **40 µL BFL1 Additive**. Vortex briefly.
- Add **20 µL Proteinase K (20 mg/ mL)**. Mix well. Incubate at **65°C, 30 minutes**. Mix 2-3 times by inversion during incubation.
- Centrifuge at **13,000 g** for **5 minutes**. Transfer the supernatant to a fresh tube.
- [Optional]* Add **10 µL RNase-A (20 µg/ µL)** and incubate for **5 minutes** at room temperature.
- Add **200 µL BFL2 Buffer** and mix by inversion.
Check isopropanol has been added to BFL2 Buffer (✓).
- Transfer **up 700 µL mixture** to a High-Q™ Spin Column placed into a Collection Tube. Incubate **2 minutes** at room temperature.
- Centrifuge at **10,000 g** for **1 minute**. Remove the flow-through and place back the High-Q™ Spin Column into a Collection Tube. Repeat steps 7 and 8 with the remaining mixture.
- Add **500 µL BFL3 Buffer**.
Check absolute ethanol has been added to BFL3 Buffer (✓).
- Centrifuge at **10,000 g** for **1 minute**. Discard the flow-through.
- Place back the High-Q™ Spin Column into the Collection Tube and add **500 µL BFL4 Buffer**. Repeat this step one more time.
Check absolute ethanol has been added to BFL4 Buffer (✓).
- To dry High-Q™ Spin Column and eliminate residual ethanol, centrifuge again at **10,000 g** for **1 minute**.
- Place the High-Q™ Spin Column into a clean 1.5 mL Tube.
- Add **50-100 µL prewarmed Elution Buffer or Water** (Molecular Biology Grade).
Prewarm Elution Buffer or Water at 70°C.
- Incubate at room temperature, **2 minutes**.
- Centrifuge at **10,000 g** for **1 minute** to elute purified DNA.
- Check DNA quantity by spectrophotometry and quality on agarose electrophoresis gel.
- Store DNA at **-20°C**.

