

# High-Q™-Spin-Column Plant Genomic DNA Purification Kit

## Ordering info

TBK0166, 3 reactions (sample)

TBK0167, 50 reactions

TBK0168, 100 reactions

TBK0169, 200 reactions

## Description

High-Q™-Spin-Column Plant Genomic DNA Purification Kit is an easy silica-membrane-based system for DNA purification from a wide variety of plant species. An optimized lysis buffer guarantees a good yield while the use of High-Q™ Spin Columns allow a good quality DNA, suitable for downstream applications.

## Features

- Starting material up to 100 mg of fresh material and up to 30 mg of dried plant material.
- Typical yields are 2- 50 µg of DNA depending on the material plant used.
- No organic extraction, no ethanol precipitation.
- High DNA purity, the isolated DNA is ready to use for downstream molecular biology applications.
- Easy and Fast protocol.

## Applications

- Purification of DNA from plant tissue, including plant cells, leaves, seeds, fruits or roots.
- Purification of DNA plant using different starting plant materials: frozen, fresh or dried.
- DNA obtained is suitable for downstream molecular biology applications such as PCR, enzymatic digestion for cloning or Southern, genotyping, etc.

## Quality Control

DNA purified is checked by: integrity (agarose gel electrophoresis), quantity and quality ( $A_{260/280} = 1.8 \pm 0.2$ ;  $A_{260/230} = 2.0 \pm 0.2$ ).

## Kit Components

| Components                                | TBK0167            | TBK0168                | TBK0169              |
|-------------------------------------------|--------------------|------------------------|----------------------|
| High-Q™ Spin Column with Collection Tubes | 50                 | 100                    | 200                  |
| BPL1 Buffer                               | 25 mL              | 47 mL                  | 92 mL                |
| BPL2 Buffer                               | 2 mL <sup>a</sup>  | 4 mL <sup>b</sup>      | 7 mL <sup>c</sup>    |
| BPL3 Buffer                               | 12 mL <sup>d</sup> | 24 mL <sup>e</sup>     | 44 mL <sup>f</sup>   |
| BPL4 Buffer                               | 12 mL <sup>g</sup> | 24 mL <sup>h</sup>     | 50 mL <sup>i</sup>   |
| Elution Buffer                            | 15 mL              | 15 mL                  | 25 mL                |
| Proteinase K                              | 30 mg <sup>j</sup> | 2 x 30 mg <sup>j</sup> | 3x30 mg <sup>j</sup> |
| Proteinase K Resuspension Buffer          | 1.5 mL             | 2 x 1.5 mL             | 3 x 1.5 mL           |

**Order Info Kit Components:** High-Q™ Spin Column with Collection Tubes (TBM0010) | BPL1 Buffer (TBB0530) | BPL1 Additive (TBB0531) | BPL2 Buffer (TBB0532) | BPL3 Buffer (TBB0533) | BPL4 Buffer (TBB0534) | Elution Buffer (TBB0510) | Proteinase K (TBZ0305) | Proteinase K Resuspension Buffer (TBB0546).

**Components for samples are ready to use!**

## Before its use:

- <sup>a</sup> Add 18 mL isopropanol and mix well.
- <sup>b</sup> Add 36 mL isopropanol and mix well.
- <sup>c</sup> Add 63 mL isopropanol and mix well
- <sup>d</sup> Add 18 mL absolute ethanol and mix well.
- <sup>e</sup> Add 36 mL absolute ethanol and mix well.
- <sup>f</sup> Add 66 mL absolute ethanol and mix well
- <sup>g</sup> Add 48 mL absolute ethanol and mix well.
- <sup>h</sup> Add 96 mL absolute ethanol and mix well.
- <sup>i</sup> Add 200 mL absolute ethanol and mix well.
- <sup>j</sup> Add 1.5 mL Proteinase K Resuspension Buffer and mix well.

## Storage

Store the kit at 25°C and Proteinase K at -20°C.

## Material required (not supplied)

- 1.5 mL Tubes.
- Ethanol (CAS 64-17-5).
- Isopropanol (CAS 67-63-0).

## PROTOCOL

1. Grind up to 100 mg of plant material in liquid nitrogen using a mortar and a pestle. With a freeze spatula, collect tube. Commercially available equipment for homogenization also can be used with zirconia beads.  
*Grind up to 50 mg if the plant material is seed.*
2. Add **450 µL BPL1 Buffer**. Vortex vigorously for 1 minute.
3. Add **20 µL Proteinase K (20 mg/mL)**. Vortex briefly. Incubate at 65°C, 30 minutes. Mix 2-3 times during incubation.
4. Centrifuge at 13,000 g for 5 minutes. Transfer the supernatant to a fresh tube.
5. **[Optional]** Add **20 µL RNase-A 10 mg/ mL** and incubate for 5 minutes at room temperature.
6. Add **300 µL BPL2 Buffer** and mix by inversion.  
*Check isopropanol has been added to BPL2 Buffer (✓).*
7. Transfer **up 700 µL mixture** to a High-Q™ Spin Column placed into a Collection Tube.
8. Centrifuge at 10,000 g for 1 minute. Remove the flow-through and place back the High-Q™ Spin Column into a Collection Tube. If there is still any remaining mixture, repeat steps 7 and 8.
9. Add **500 µL BPL3 Buffer**.  
*Check absolute ethanol has been added to BPL3 Buffer (✓).*
10. Centrifuge at 10,000 g for 1 minute. Discard the flow-through and place back the High-Q™ Spin Column into a Collection Tube.
11. Add **500 µL BPL4 Buffer**. Repeat this step one more time.  
*Check absolute ethanol has been added to BPL4 Buffer (✓).*
12. To dry the High-Q™ Spin Column and eliminate residual ethanol, centrifuge again at 10,000 g for 2 minutes.
13. Place the High-Q™ Spin Column into a clean 1.5 mL Tube.
14. Add **50-100 µL prewarmed Elution Buffer or Water (Molecular Biology Grade)**.  
*Prewarm Elution Buffer or Water at 70°C.*
15. Incubate at room temperature, 2 minutes.
16. Centrifuge at 10,000 g for 1 minute to elute purified DNA.
17. Check DNA quantity by spectrophotometry and quality on agarose electrophoresis gel.
18. Store DNA at -20°C.