

STOUT™ FAST PCR Master Mix (2x)

Ordering Info

TBK0051, 80 reactions

TBK0052, 400 reactions

Description

STOUT™ FAST PCR Master Mix (2x) is a convenient formulation to daily fast PCR reactions. It includes STOUT™ FAST DNA polymerase a recombinant enzyme with a fast polymerization range (4-8 kb/ min) generating consistent amplification results. The enzyme has 5'→ 3' polymerase activity and 5'→ 3' exonuclease activity while 3'→ 5' exonuclease activity is absent. It is an extraordinary enzyme for routine PCR of genotyping or screening.

Features

- **Fast Amplification** of PCR targets up to 5 kb.
- **High Extension Rate**, 2 seconds/ kb for targets < 1 kb.
- **Addition without template** of 3'adenine at the end of PCR fragment.
- **Non-proofreading** polymerase.
- **Immediate Activation**.

Applications

- Fast PCR.
- Generation of PCR fragments for TA cloning.
- Genotyping.
- Screening by PCR.

Kit Components

Components	TBK0051	TBK0052
STOUT™ FAST PCR Master Mix (2x)	1 mL	5 x 1 mL
PCR Grade Water	2 mL	3 x 2 mL

Order Info Kit Components: STOUT™ FAST Master Mix (TBK0051) | PCR Grade Water, nuclease free (TBK0303).

Storage

Store at -20°C. Shipped in blue ice.

Quality Control

Functionally tested in a 1 kb PCR amplification (GC 52%).

Material required (not supplied)

- PCR Tubes
- Specific Primers

Also available:

- **STOUT™ FAST DNA Polymerase** (TBK1011, TBK1012)
- **STOUT™ FAST RED Master Mix 2x** (TBK0053, TBK0054)
- **WARM™ FAST Hot-Start Red Master Mix (2x)** (TBK0039, TBK0040)
- **STOUT™ Recombinant Taq DNA Polymerase Master Mix (2x)** (TBK0028, TBK0029)

PROTOCOL

1. Thawing all components on ice. Vortex and centrifugate them.
2. On ice, prepare a mix of the following components, considering the number of samples plus two extra reactions.

Reaction Components	Final Concentration	Volume
STOUT™ FAST PCR Master Mix 2x	1 x	12.5 µL
Forward Primer (10 pmol/ µL)	0.4 µM	1 µL
Reverse Primer (10 pmol/ µL)	0.4 µM	1 µL
PCR Grade Water		up 25 µL*
DNA template (add in step 4)		*
Final Volume		25 µL

* Consider volume of template to be added in step 4.

3. Distribute the mix prepared in each PCR tube or well.
4. Add in each tube the DNA sample (cDNA: < 50 ng; gDNA: < 250 ng). Mix well.
5. Set up thermocycler:

Process	Cycles	Temperature	Time
Initial denaturation	1 x	94 °C	1:00
Denaturation	25-40 x	94 °C	0:15
Annealing		Tm	0:15
Extension		72 °C	0:02 per kb (< 1 kb) 0:15 per kb (> 1 kb)
Final Extension	1 x	72 °C	3:00
Conservation	1 x	4 °C	∞