

Senescence Detection Kit

(Colorimetric Readout)

Ordering info

TBK0525, 100 assays

Description

Senescence Detection Kit enables reliable identification of senescent cells *in vitro* by detecting the enzymatic activity of lysosomal senescence-associated β -galactosidase (SA- β -gal), a widely accepted specific biomarker of cellular senescence. This assay is based on the hydrolysis of substrate (5-bromo-4-chloro-3-indolyl- β -D-galactopyranoside), a chromogenic substrate, which yields an easily visible blue precipitate under light microscopy.

β -Galactosidase activity is significantly increased in senescent cells, allowing their clear identification even in heterogeneous cultures. The enzyme hydrolyzes the substrate into a colored product that accumulates within the cells, enabling visual detection of them.

Features

- **Direct visual detection of senescence**, through the formation of a blue precipitate.
- **Compatible** with post-treatment analysis.
- **Highly specific**, differentiation between SA- β -gal (pH 6.0) and basal lysosomal activity.

Applications

- Evaluation of anti-senescence therapies (senolytics).
- Characterization of cells at the Hayflick limit.
- Phenotypic profiling of tumor or primary cell cultures.

Kit Components

Components	TBK0525
Substrate X-Gal ^a	6 x 20 mg
Fixative Solution 10x ^b	15 mL
Staining Solution 5x ^c	15 mL
Additive Solution 1, 100x	1.5 mL
Additive Solution 2, 100x	1.5 mL
PBS 1x, pH 7.4 (powder) ^d	1 pouch

Order Info Kit Components: Substrate (TBK0525-1) | Fixative Solution 10x Buffer (TBK0525-2) | Staining Solution 5x (TBK0525-3) | Additive Solution 1 (TBK0525-4) | Additive Solution 2 (TBK0525-5) | PBS pouch (TBK0600-1)

Before its use

^a **Substrate X-Gal Solution 20x:** Dissolve **20 mg Substrate** in 1 mL N,N, Dimethylformamide (CAS 68-12-2). Use a glass or a polypropylene container for X-GAL. Store at -20 °C.

^b **Fixative Solution 1x:** Dilute 1 part of **Fixative Solution 10x** in **9 parts of Water**. Store at -20 °C.

Fixative solution is toxic; it contains formaldehyde and glutaraldehyde and should therefore be handled with care to avoid inhalation or skin exposure.

^c **Staining Solution 1x:** Dilute 2 part of **Staining Solution 5x** in 8 parts of Water. Store at 4°C.

^d **PBS 1x:** Dissolve 1 pouch of PBS pH 7.4 in 1 L of sterile distilled water. Store at room temperature.

Staining Working Solution pH 6.0: In a polypropylene tube add,

	To prepare	
	1 mL	10 mL
Staining Solution 1x	930 μ L	9.3 mL
Additive Solution 1, 100x	10 μ L	100 μ L
Additive Solution 2, 100x	10 μ L	100 μ L
Substrate Solution 20x	50 μ L	500 μ L

Storage

Store the kit protected from light at -20 °C.

PROTOCOL

This protocol is optimized for individual wells of a 6-well plate and can be adjusted as needed for other culture plate formats.

1. Seed $2-5 \times 10^4$ cells in either a 35-mm plate or 6-well plate, and culture for 2–3 days or more.
2. Remove culture medium and wash cells with **2 mL PBS 1x**.
3. For fixing cells or frozen tissue sections, add **1-2 mL Fixative Solution 1x** to submerge the cells.
4. Incubate for 10-15 minutes at room temperature.
5. Remove **Fixative Solution 1x** and wash the cells twice with **2 mL PBS 1x**.

6. Add **1 mL Staining Working Solution pH 6.0** to each well and cover the plate to prevent evaporation.

The pH of the Staining Working Solution must be 6.0 to accurately detect SA- β -gal activity. A lower pH may lead to false positives, while a higher pH may cause false negatives. If necessary, adjust the pH using HCl or NaOH.

7. Incubate in the dark at 37°C, for 2 hours up to overnight.

Do not incubate in a CO₂ incubator. CO₂ can lower the pH and interfere with the development of the colorimetric reaction. Blue color is detectable in some cell lines within 2 hours, but staining is maximal at 12-16 hours.

8. After blue color is fully developed wash cells twice with **PBS 1x**.
9. Add one drop of mounting medium, and place cover glasses either on a 35-mm plate or 6-well plate. Observe the cells under a microscope.
10. Count the number of blue stained senescent cells and total cells. Calculate the percentage of senescent cells.