

WARM™ Fast Hot-Start DNA Polymerase Kit

Ordering Info

TBK1026, 250 U

TBK1027, 500 U

Description

WARM™ Fast Hot-Start DNA Polymerase Kit is a convenient kit to speed amplification of DNA samples. The kit includes the robust WARM™ Fast Hot-Start polymerase, characterized by its extension rate of 2 sec/ kb and an optimized buffer composed of MgCl₂ 15 mM, highly purity High-Q™ dNTPs 5 mM each, and enhancers. WARM™ Hot-Start DNA polymerase is a modified polymerase suitable for amplification from unpurified samples or difficult targets.

Features

- **Amplification up 5 kb** amplicons.
- **Extension rate** 2 seconds per kb (< 1kb) and 15 seconds per kb (> 1 kb to 5 kb).
- **Suitable for TA cloning purposes**, the enzyme adds 3' adenine at the end of PCR fragment.
- **Activation Controlled**, inactive at low temperatures and fully activated at high temperature.

Applications

- Amplification of complex targets.
- Direct PCR of unpurified samples (*blood, urine*).
- Generation of PCR fragments for TA cloning.
- Multiplex PCR.

Kit Components

| Components | TBK1026 | TBK1027 |
|--|----------|----------|
| WARM™ Fast Hot-Start DNA Polymerase (5 U/μL) | 50 μL | 2 x 50μL |
| 5x WARM™ Fast Hot-Start Buffer | 2 x 1 mL | 4 x 1 mL |
| PCR Grade Water, nuclease free | 2 mL | 2 x 2 mL |

Order Info Kit Components: WARM™ Fast Hot-Start DNA Polymerase (TBK1026-1) | 5x WARM™ Fast Hot-Start Buffer (TBK1026-2) | PCR Grade Water, nuclease free (TBB0303).

Storage

Store at -20°C protected from light. Shipped in blue ice.

Quality Control

Functionally tested in a PCR amplification.

Material required (not supplied)

- PCR Tubes
- Primers

PROTOCOL

I. PREPARING REACTIONS

1. Gently vortex and briefly centrifuge all solutions after thawing.
2. Place a tube on ice and add the following components for each 25 μL reaction. Prepare sufficient master mix for the number of reactions. Consider one or two extras:

| Components | Volume | Final Concentration |
|---|------------------------------------|---------------------|
| 5x WARM™ Fast Hot-Start Buffer | 5 μL | 1x |
| Forward primer 5 μM (5 pmol/ μL) | 2 μL | 0.4 μM |
| Reverse primer 5 μM (5 pmol/ μL) | 2 μL | 0.4 μM |
| PCR Grade Water, nuclease-free | up 25 μL | |
| Template DNA | | (1) |
| WARM™ Fast Hot-Start DNA Polymerase (5 U/ μL) | 0.1 μL | 0.5 U (2) |
| Final Volume | 25 μL | |

(1) Optimal amounts of template: genomic DNA ≤ 250 ng, cDNA synthesis reaction 1-5 μL (≤ 50 ng), direct PCR of unpurified samples from urine or blood 1 μL .

(2) Higher enzyme concentration may be required for difficult amplifications (up 2.5 U)

3. Aliquot the master mix into individual PCR tubes or in well-plate.
4. Gently vortex and spin down the samples. Add template DNA.

II. PCR SETUP

1. Perform PCR using standard thermal cycling conditions:

| Step | Cycles | Temperature | Time |
|-----------------------------------|--------|----------------|-------------|
| Activation & Initial Denaturation | 1 x | 95 °C | 3:00 |
| Denaturation | | 95 °C | 0:15 |
| Annealing | 40 x | T _m | 0:15 |
| Extension | | 72 °C | sec per kb* |
| Final Extension | 1 x | 72 °C | 3:00 |
| Conservation | 1 x | 4 °C | ∞ |

* < 1kb, 2 seconds per kb; >1-5 kb, 15 seconds per kb; multiplex PCR, 90 seconds per kb independent of size.

2. Store the samples at -20°C until use.