# TIARISTAIN™ Green Safe

## (Molecular Biology Grade)

### Ordering info

TBR0225, TIARISTAIN™ Green Safe, 20,000x, 20 µL

TBR0226, TIARISTAIN™ Green Safe, 20,000x, 1 mL

### Description

TIARISTAIN™ Green Safe is a new formulation to stain nucleic acids in agarose and polyacrylamide gels. It emits a green fluorescence when bound to nucleic acids. TIARISTAIN™ Green Safe has a strong excitation peak at 485 nm enabling the use of less DNA destructive blue light instead UV light (better low-pressure lamp, 254 nm). The dye has also two excitation peaks in the ultraviolet range, at 270 nm and 290 nm). The maximum emission is at 525 nm.

### Features

- Non carcinogenic high-quality product.
- Excited with both blue or ultraviolet light and has a maximum emission at 525 nm (green).
- High sensitivity, 0.2-0.6 ng DNA per band.
- Flexible use, valid for pre-casting and post-staining use.

### Applications

- dsDNA, ssDNA and RNA staining in agarose gels.
- Non-carcinogenic alternative to ethidium bromide nucleic acid stain.

### Storage

Store at room temperature. Protect from light.

Do not store at 2°C to 8°C. Do not freeze.

### Quality Control

DNase/ RNase activity not detected.

### Recommended devices

Transilluminators with blue light (470-520 nm), ultraviolet light (preferable with low pressure lamp 254 nm; high pressure 365 nm but excitation is less efficient).

Also available:

**INTENSE™ INDIGO™ Loading Buffer 6x (TBB0325)**

**COBALT™ Loading Buffer 6x (TBB0321)**

**GREENY™ Loading Buffer 6x (TBB0322)**

**CORAL™ Loading Buffer 6x (TBB0323)**

**COLLAGE™ Loading Buffer 6x (TBB0324)**

## PROTOCOL

; Wear gloves while handling!

### I. PRE-CASTING

1. Prepare an agarose gel (0.5 - 3% in TAE or TBE Buffer 1x).
2. Heat the solution until all the agarose is completely melted.
3. Add TIARISTAIN™ Green Safe and mix thoroughly,

V (Agarose Gel)	V (TIARISTAIN™ Green Safe, 20,000x)
100 mL	5 µL
50 mL	3 µL

4. Cool the solution and pour into the tray.
5. Once the gel is solid, load the samples.
6. For best results, add 1 µL /10 mL running buffer near the anode.

### II. POST STAINING

*Post-staining gives a better sensitivity than precast gels.*

1. Prepare the Post-Staining Solution adding **20 µL TIARISTAIN™ Green Safe 20,000x** to 100 mL TAE or TBE Buffer 1x (agarose gel thickness <0.5 cm). Mix well.

*Post-Staining Solution could be re-used 2-3 times. Store in the dark at room temperature.*

2. Soak the gel at room temperature for 10-60 minutes.

*Incubation time depends on agarose percentage and gel thickness.*

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