

ONPG β -Galactosidase Assay Kit

(Colorimetric Readout)


Ordering info

TBK0543, 500 assays

Description

ONPG β -Galactosidase Assay Kit is an extraordinary tool commonly used in molecular biology and microbiology research to measure the activity of beta-galactosidase enzyme. This assay can be used to detect the presence of certain bacterial strains and to monitor expression of *lacZ* as β galactosidase reporter gene.

ONPG β -Galactosidase Assay Kit is a quantitative method based on the cleavage of the synthetic chromogenic lactose analogous substrate ortho-nitrophenyl-beta-galactopyranoside (ONPG), releasing a bright yellow product called ortho-nitrophenol (ONP). The ONP production per unit time is directly proportional to the activity of β -galactosidase in the sample.



Features

- **Accurate, simple and reproducible** assay, around 3,5-6% inter- and 3,3-5,5% intra- assay coefficient of variation.
- **Narrow dynamic range**, linearity from 0.2 to 2 nmol/min/mg protein.
- **Colorimetric readout** at 420 nm.
- **Highly Sensitive**, the assay can detect beta-galactosidase activity in the range of 0.1-10 nmol/min/mg protein.
- **Proteolytic resistance** of β -galactosidase, lysates could be assayed directly or store at -80 °C for at least 2 months.

For Research Use Only

© TIARIS Biosciences, 2026

Kit Components

Components	TBK0543
BGA Lysis Solution ^a	50 mL
10x BGA Reaction Buffer ^a	12 mL
BGA Stop Solution ^a	70 mL
BME, β mercaptoethanol ^a	500 μ L
β galactosidase (0.4 U/ μ L) ^b	100 μ L
ONPG Solution (4 mg/ mL) ^b	10 mL

Order Info Kit Components: BGA Lysis Solution (TBB0570) | 10x BGA Reaction Buffer (TBB0571) | BGA Stop Solution (TBB0572) | Beta galactosidase (TBZ0330) | ONPG Solution (TBR0267) | BME (TBR0107).

Storage

^a Store at 4 °C

^b Store at -20 °C. Protect ONPG Solution from light.

The kit is shipped with blue ice.

Unit Definition

One unit of β -Galactosidase hydrolyzes 1 micromole of ONPG to ONP and galactose per minute at pH 7.5 and 37 °C.

Quality Control

The kit is tested in a functional assay.

Applications

- Normalization of transfection efficiency in mammalian cells.
- Study of protein trafficking, strength of promoters and enhancer with β -galactosidase as reporter.
- Protein-protein interaction studies.

PROTOCOL

I. REAGENT PREPARATION

- **2x Reaction Buffer:** Prepare enough buffer to your assay setup. For 100 reactions, dilute 1 mL 10x BGA Reaction Buffer with 4 mL Deionized Water. Add 17,5 μL β -mercaptoethanol. This buffer could be stored at 4 °C.

II. CELL LYSATE PREPARATION

For transient expression, cell lysates must be prepared 48-72 hours after plasmid transfection. As negative control, prepare cell lysates from the same cell line that have not been transfected.

Suspension cells

1. Centrifuge the cells (~ 2 mL) at 250 g, 5 minutes. Eliminate growth medium.
2. Wash twice the cells with PBS 1x (without Ca^{2+} and Mg^{2+}). Discard the supernatant.
3. Add 50 μL BGA Lysis Solution and incubate at 37 °C, 30 minutes.

Adherent cells

1. Discard growth medium and wash once the cells with PBS 1x (without Ca^{2+} and Mg^{2+}). Discard PBS 1x.
2. Add a sufficient volume of BGA Lysis Solution to cover the cells. The next table could be indicative.

Culture Format	BGA Lysis Solution per well
96-well plate	50 μL
24-well plate	250 μL
12-well plate	500 μL
6-well plate	800 μL
60 mm dish	1500 μL

If transfection efficiency is very high, it may be necessary to dilute the lysate in BGA Lysis Solution. On the contrary, if the expression is low, reduce the volume of BGA Lysis Solution.

3. Incubate at room temperature for 15 minutes, swirling the plates several times to ensure an optimal lysis.
Optimal lysis is essential to obtain reproducible results.
4. Transfer the cell lysate to a microcentrifuge tube by pipetting.
5. Spin the lysates at 4 °C for 2 minutes to eliminate cellular debris.
6. Transfer the supernatant to a microcentrifuge tube.

The lysates may be assayed directly or stored at -80 °C for at least 2 months.

III. ASSAY SETUP

All samples must run at least in duplicate, preferably triplicate. Keep the samples and standard on ice. Equilibrate at room temperature the rest of components.

- **Standard curve:** Firstly, prepare fresh β galactosidase Solution at 40 mU/ μ L. Make 1:10 dilution of the provided enzyme, adding 10 μ L β -galactosidase (0.4 U/ μ L) to 90 μ L BGA Lysis Solution. After that, in clean tubes, prepare a standard curve as follow:

β -galactosidase Activity (U/ μ L)	10	8	6	4	2	1	0.5	0
β -galactosidase Solution 40 mU/ μ L (μ L)	10	8	6	4	2	1	0.5	0
BGA Lysis Solution (μ L)	30	32	34	36	38	39	39.5	40

Add 10 μ L of each dilution to its corresponding wells of 96-well plate.

- **Negative Control:** Cells lacking β -galactosidase, such as LT2 or non-transfected cells.
- **Positive Control:** Cells with high β -galactosidase activity.
- **Plate setup:**

	1	2	3	4	5	6	7	8	9	10	11	12
A	10 U/ μ L			Sample 1			Sample 8			Sample 15		
B	8 U/ μ L			Sample 2			Sample 9			Sample 16		
C	6 U/ μ L			Sample 3			Sample 10			Sample 17		
D	4 U/ μ L			Sample 4			Sample 11			Sample 18		
E	2 U/ μ L			Sample 5			Sample 12			Sample 19		
F	1 U/ μ L			Sample 6			Sample 13			Sample 20		
G	0.5 U/ μ L			Sample 7			Sample 14			Sample 21		
H	0 U/ μ L			NEGATIVE CONTROL			POSITIVE CONTROL			Sample 22		

STANDARD CURVE
 BLANK
 NEGATIVE CONTROL
 POSITIVE CONTROL
 SAMPLES

IV. BETA GALACTOSIDASE ASSAY

- Following the plate setup:
 - Add 10 μ L of dilution prepared for **Standard Curve**.
 - Add 10 μ L **Cell Lysate** from samples and from negative and positive control.
- With a multichannel pipette, add 15 μ L **ONPG Solution** to each well.
- With the multichannel pipette, add 50 μ L **1x BGA Reaction Buffer** to each well.
- Shake the plate carefully for 30 seconds.
- Cover the plate with metallic film and incubate it at 37 °C for 30 minutes.
- Add 125 μ L **BGA Stop Solution**
- Read the absorbance at 420 nm.

V. BETA GALACTOSIDASE ACTIVITY DETERMINATION

The absorbance from the blank control well is used as the background absorbance. Subtract the background absorbance from samples absorbance.

To quantify the amount of β -galactosidase in the samples, interpolate them on the enzyme standard curve.

For Research Use Only

© TIARIS Biosciences, 2026

ONPG β -Galactosidase Assay Kit