

TIARIS™ Ruler 100 – 2.000 bp

(Molecular Biology)


Ordering info

TBR0239, 50 µg, 0.1 µg/ µL

TBR0240, 5 x 50 µg, 0.1 µg/ µL

Description

TIARIS™ Ruler 100-2000 bp is a ready to use molecular weight marker. It is composed of 13 fragments ranging from 100 bp to 2000 bp. 500 bp band has increased intensity relative to the other bands for easy identification. The product includes loading buffer with a loading dye.



Features

- Ready to use.
- Bright and accurate.
- Included 500 bp-highlighted band as visual reference.
- Suitable for TAE and TBE Buffer.

Storage

Store at -20°C.

The product is shipped at room temperature (stable 3 months).

Applications

- Size determination of restriction digests, PCR products, DNA fingerprinting profiles, etc.
- Quantitation of dsDNA.

Quality Control

Band location and intensity by agarose gel electrophoresis.

Also available:

STOUT™ RED PCR Master Mix (2x) (TBK0026)

TAE 10x, pH=8.3 (TBB0355)